


High Performance Liquid Chromatography (HPLC) - Purity Test



Salts and Sugars used as buffering agents are not treated as impurities


Chromatogram Peaks

Number of detected peaks: 5

Peak.	RT.	Item.
1.	1.740	GHK-CU
2.	8.360	BPC-157
3.	2.097	TB-500
4.	1.464	KPV
5.	3.814	

Analytical methods: UV (214-290 nm), HPLC (Purity), and LC-MS (Identity)
Combined Purity: 99.96%
Compound: GHK-CU, BPC-157, TB-500, Msh (11-13)
Pubchem CID: 71587328, 9941957, 62707662, 125672
Analysis date: 03/02/26
Calculated Molecular Mass - 402.9, 1419.5, 889.0, 342.4 Da
Quantity: 51.31, 10.84, 10.17, 10.05 mg

Mass Spectrometric Analysis



Analysis performed by:
Richard Sylvester MD

DR Richard Sylvester

LOT-KLOW-00003