





High Performance Liquid Chromatography (HPLC) - Purity Test




Salts and Sugars used as buffering agents are not treated as impurities

Chromatogram Peaks			Number of detected peaks: 4
Peak.	RT.	Item.	Analytical methods: UV (214-290 nm), HPLC (Purity), and LC-MS (Identity) Purity: 99.23% Compound: Pramipexole Pubchem CID: 119570 Analysis date: 03/21/26 Calculated Molecular Mass - 211.3 Da
1.	2.694	Pramipexole	
2.	2.189		
3.	3.781		
4.	4.796		

Mass Spectrometric Analysis



Analysis performed by:
 Richard Sylvester MD



LOT-PRAM-00001