




High Performance Liquid Chromatography (HPLC) - Purity Test



Salts and Sugars used as buffering agents are not treated as impurities


Chromatogram Peaks

Number of detected peaks: 2


Peak.	RT.	Item.
1.	11.661	Sildenafil
2.	2.189	

Analytical methods: UV (214-290 nm), HPLC (Purity), and LC-MS (Identity)
Purity: 99.12%
Compound: Sildenafil
Pubchem CID: 135398744
Analysis date: 03/02/26
Calculated Molecular Mass - 474.6 Da

Mass Spectrometric Analysis



Analysis performed by:
Richard Sylvester MD



LOT-SILD-00002