




High Performance Liquid Chromatography (HPLC) - Purity Test



Salts and Sugars used as buffering agents are not treated as impurities


Chromatogram Peaks

Number of detected peaks: 2

Peak.	RT.	Item.
1.	12.332	SLU-PP-332
2.	5.685	

Analytical methods: UV (214-290 nm), HPLC (Purity), and LC-MS (Identity)
 Purity: 99.91
 Compound: SLU-PP-332
 Pubchem CID: 5338394
 Analysis date: 11/15/25
 Calculated Molecular Mass - 290.2 Da

Mass Spectrometric Analysis



Analysis performed by:
 Richard Sylvester MD



LOT-SLUP-00009